



Letter code	Tube Size in mm	13/75				13/100		16/100		9/120	
		2	3	3.5	4	5	6	8	9	1.5	2.75
Z	Serum Clot Activator	PREM/NR	PREM/NR		PREM/NR		PREM/NR		NR		
	Serum Clot Activator + Gel				PREM/NR	PREM/NR		NR	NR		
LH	Lithium Heparin	PREM/NR	PREM/NR		PREM/NR		PREM/NR		NR		
	Lithium Heparin + Gel		NR		PREM/NR	PREM/NR		NR			
NH	Sodium Heparin				PREM/NR		PREM/NR		NR		
9NC	Trisodium citrate solution	PREM/NR	PREM/NR	PREM/NR							
K3E	K3 EDTA	PREM/NR	PREM/NR	PREM/NR	PREM/NR	PREM	PREM/NR		NR		
K2E	K2 EDTA	PREM/NR	PREM/NR		PREM/NR	PREM/NR	PREM/NR		NR		
K2E	K2 EDTA + Gel				PREM/NR	PREM/NR		NR			
FX	Sodium fluoride Potassium oxalate		PREM		PREM/NR		PREM				
FE	Sodium fluoride EDTA K3	PREM/NR			PREM/NR						
ESR	Trisodium citrate	PREM									PREM

PREM = PREMIUM (Safety Twist Cap), NR = Non-ridged (pull cap)

**VACUETTE® Coagulation**

These tubes are used for determinations in citrated plasma for coagulation testing. The recommended, buffered sodium citrate solution 3.2% (0.109mol/l) functions as anticoagulant by chelating calcium. The proportion of blood to sodium citrate anticoagulant volume is 9:1.

**VACUETTE® Serum**

Serum Tubes are used for determinations in serum for clinical chemistry, microbiological, serology, immunology and TDM. The tubes are coated with a clot activator. The recommended clotting time is >=30 minutes. Serum Gel Tubes allow storage of certain parameters under the recommended storage conditions according to their biological half-life in separated serum.

**VACUETTE® Heparin Plasma**

The tubes are used for determinations in heparinised plasma for clinical chemistry. The heparin concentration is standardised: 18 I.U. of lithium-, sodium- or ammonium salt of heparin per 1ml. Plasma Gel Tubes allow storage of certain parameters under the recommended storage conditions according to their biological half-life in separated plasma.

**VACUETTE® EDTA**

EDTA Tubes are used for determinations in EDTA whole blood for haematology and immunohaematology. The tube interior of EDTA Tubes is spray dried with 1.2 - 2mg anhydrous EDTA per 1ml blood.

**VACUETTE® Glucose**

Glucose Tubes are used for the analysis of blood sugar and lactate in stabilised whole blood or plasma. Sodium fluoride is used as glycolysis inhibitor to preserve glucose when combined with an anticoagulant such as potassium oxalate or EDTA.

**VACUETTE® ESR**

Sodium citrate solution 3.2% (0.109 mol/l) is the anticoagulant of choice for the measurement of the ESR. The proportion of blood to sodium citrate anticoagulant volume is 4:1.



For further information, please visit our website [www.gbo.com/preanalytics](http://www.gbo.com/preanalytics) or contact us:

**Austria (Main Office)**  
Greiner Bio-One GmbH  
Phone +43 7583 6791-0  
Fax +43 7583 6318  
E-Mail [office@at.gbo.com](mailto:office@at.gbo.com)

**Germany**  
Greiner Bio-One GmbH/Preanalytics  
Tel +49 7022 948-0  
Fax +49 7022 948-514  
E-Mail [office@de.gbo.com](mailto:office@de.gbo.com)

**Spain**  
VACUETTE Espana S.A.  
Phone +34 91 652 77 07  
Fax +34 91 652 33 35  
E-Mail [info@vacuette.es](mailto:info@vacuette.es)

**United Kingdom**  
Greiner Bio-One Ltd.  
Phone +44 1453 8252 55  
Fax +44 1453 8262 66  
E-Mail [info@uk.gbo.com](mailto:info@uk.gbo.com)

**Brazil**  
Greiner Bio-One Brasil  
Phone +55 19 3468-9600  
Fax +55 19 3468-9621  
E-Mail [office@br.gbo.com](mailto:office@br.gbo.com)

**Hungary**  
Greiner Bio-One Hungary Kft.  
Phone +36 96 213 088  
Fax +36 96 213 198  
E-Mail [office@hu.gbo.com](mailto:office@hu.gbo.com)

**Switzerland**  
Greiner Bio-One VACUETTE Schweiz GmbH  
Phone +41 7 12 28 55 22  
Fax +41 7 12 28 55 21  
E-Mail [office@ch.gbo.com](mailto:office@ch.gbo.com)

**USA**  
Greiner Bio-One North America Inc.  
Phone +1 704 261-7800  
Fax +1 704 261-7899  
E-Mail [info@us.gbo.com](mailto:info@us.gbo.com)

**China**  
Greiner Bio-One Suns Co., Ltd.  
Phone +86 10 83 55 19 91  
Fax +86 10 63 56 69 00  
E-Mail [office@cn.gbo.com](mailto:office@cn.gbo.com)

**India**  
Greiner Bio-One INDIA Pvt., Ltd.  
Phone +91 120 456 8787  
Fax +91 120 456 8788  
E-Mail [info@gboindia.com](mailto:info@gboindia.com)

**Thailand**  
Greiner Bio-One Thailand Ltd  
Phone +66 38 4656 33  
Fax +66 38 4656 36  
E-Mail [office@th.gbo.com](mailto:office@th.gbo.com)

**France**  
Greiner Bio-One SAS  
Phone +33 1 69 86 25 25  
Fax +33 1 69 86 25 35  
E-Mail [office@fr.gbo.com](mailto:office@fr.gbo.com)

**Netherlands**  
Greiner Bio-One B.V.  
Phone +31 172 4209 00  
Fax +31 172 4438 01  
E-Mail [info@nl.gbo.com](mailto:info@nl.gbo.com)



980044 rev.02, 06.2014 e



**VACUETTE® Specimen Collection System**

Our Innovations for Your Safety



Thread the needle into the needle holder. Ensure that the needle is firmly seated in the holder. Leave the remaining section of the cover in place.



Apply tourniquet if required and disinfect the venipuncture site.



Place the patient's arm in a downward position and insert the needle into the vein as usual.



Push the tube into the holder and onto the needle valve puncturing the rubber stopper. Hold in place with the thumb. The prespecified vacuum of the tube allows the required quantity of blood to flow into the tube.



Remove the tourniquet as soon as blood appears in the tube. When the first tube is full and the blood flow ceases, remove it from the holder. Further samples can then be collected by inserting more tubes into the holder.



In order to ensure correct mixing of the sample with the additives, gently invert the tubes 5-10 times immediately following collection. After the last tube has been drawn and the blood flow ceases, remove the needle from the vein.



Activate the safety mechanism with the aid of solid surface or thumb upon removal of the needle from the patient's vein. An audible click indicates that the safety shield has been fully activated.



Dispose of the needle together with holder in a sharps disposal container specially intended for this purpose.

### Components


Multiple Use Drawing Needle or VISIO PLUS Needle



VACUETTE QUICKSHIELD or Standard Tube Holder



VACUETTE Tourniquet or VACUETTE Super-T



VACUETTE Blood Collection Tube



Disposal Containers


## VACUETTE® Specimen Collection System



The invention of the world's first evacuated blood collection system made out of PET plastic by Greiner Bio-One has made specimen collection a lot safer. The innovative VACUETTE® system guarantees simple handling and hygiene in order to make your daily work easier.

for more details please see:  
[www.gbo.com/preanalytics](http://www.gbo.com/preanalytics)

### Recommended order of draw:



\* Coagulation tubes may be the first tube to be drawn for routine testing only (PT und aPTT).  
Note: Always follow your facility's protocol for order of draw.

The VACUETTE® system combines the advantages of vacuum technology with unique safety characteristics for the patient and the user.

### Specification

- clear as glass but **virtually unbreakable**
- **screw caps** for **increased safety** and **simplified opening**
- **colour coded** rings offer additional **visual identification** of tube characteristics
- **longer shelf-life** due to **tube design**
- **compatible** with **all common analyser systems**
- **Fill-line** on each tube for **simple visual check**

### Safety

- Virtually **unbreakable** during centrifugation, handling and transportation
- all **VACUETTE® tubes are sterile**
- **reproducibility of test results is guaranteed**
- **fast & reliable test results**
- **VACUETTE® Sandwich tube** with **special double wall technology**; the ultimate standard for coagulation tests
- **incineration** of **VACUETTE® tubes is environmental friendly** as PET breaks down into 3 natural components: carbon, hydrogen and oxygen

### Standards

- **VACUETTE® products are US FDA approved**
- manufacturing complies with EC directives and US regulations (ISO 9001, ISO 13485, GMP)
- Products can be used according to recommendations by WHO and CLSI