



Item number 1000021314

ABENA®

ABENA San, Special, light green colourcode, Premium

- ✓ Designed for faecal incontinence
- ✓ High performance core with Top-Dry Technology
- ✓ Unique leakage barriers with built-in pockets
- ✓ Support healthy skin balance
- ✓ Odour Control



Product Description

ABENA San Special is our anatomically shaped pad specifically designed for faecal incontinence. The pad can be used for both faecal and urinary incontinence and features a light green textile-like breathable backsheet foil for optimum comfort and discretion. The pad has effective leakage barriers with built-in pockets that keep the faecal contents inside.

The high performance core with Top-Dry Technology ensures fast absorption and a dry surface even after multiple voids, maintaining skin integrity. The fully breathable pads are fitted with standing leakage barriers combined with front and rear built-in cross barriers offering 360° leakage protection. When combined, these functionalities keep the skin in balance and maintains a natural pH skin balance.

The odour control system minimises unpleasant smells. The wetness indicator guides the user and caretaker to the product with the right absorption level and product changes. The dermatologically tested products feature the Nordic Eco Label.

ABENA San Special should always be used with special fixation underwear, like our Abri-Fix products, to ensure a close fit and reduce the risk of leakage.



Specifications

Brand	ABENA
Sub-brand	San
Size	Special
Color	White
Features	Light green colourcode
Material	Breathable textile backsheet (BTBS), ECF fluff, elastane, nonwoven, PE, PP, resin, SAP
Length/depth	730 mm
Width	370 mm
Tolerance	+/- 10%
Absorption Rothwell	2000 ml
Gender	Unisex
Quality level	Premium
Certifications	CE, MD, Nordic Swan Ecolabel, Dermatologically tested, UKCA
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1, ISO 13485, ISO 14001, ISO 9001
Directives, regulations and acts	MDR (EU) 2017/745
Safety Instructions and Warnings	Keep out of reach of children.
Because of the products composition it has an expiry date from the production date on:	5 years

Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Dispose with household waste.
Packaging Disposal Instructions	Can be recycled or incinerated.

Instructions for use/application

Double incontinence pad designed for combined faecal and urinary incontinence. Easy product recognition by coloured backsheet foil including printed naming.

ABENA San Special should always be used with special fixation underwear, like our Abri-Fix products. The fixation pants ensures a close fit and reduces the risk of leakages. Pull fixation pants half-way up the legs. Fold the pad lengthways to create a bowl shape. Pass the pad between the legs from front to back. Position the upper rear end of the pad at the coccyx. Fan out the back of the pad horizontally to create a channel shape between the legs. Use the wetness indicator to check if the product is applied accurately in the center. Fan out the pad horizontally at the front. Pull fixation pants up into position. See our application guide brochure for detailed information.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cII	120 pcs	60 cm	40 cm	45.2 cm	5713571000717
pck	30 pcs	34.5 cm	20 cm	29 cm	5713571000700



Nordic Swan Ecolabel is an official Scandinavian environmental label. Products with the Nordic Swan Ecolabel meet stringent environmental requirements with regard to raw material, production, use, documentation, and safety.



Your guarantee that our products have been tested for ingredients which can potentially cause skin irritation. Tests conducted by an independent third-party laboratory.



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



The product does not contain latex.



Leakage barriers in the leg cut and around the absorbent core help prevent leakages and provide an optimal fit.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.



The unique structure of our high performance core is designed for fast dispersion and encapsulation of urine, which results in a dry surface even after multiple voids.



Odour control is a critical factor in ensuring discreteness for the user. ABENA's incontinence products are designed and engineered to ensure that minimal odour arises from the product.



Unique Top-Dry Technology that provides a rapid absorption and ensures a dry surface even after multiple voids.



Wetness indicator is a liquid-sensitive strip that changes colour or fades out indicating when it is time to change the product.



A combination of 100% breathability, Top-Dry and leakage protection provides a dry surface and healthy skin even after multiple voids maintaining the skin's natural pH balance.



Drop system indicating absorption. The more drops, the higher the absorption. Ranges from 1 to 9 drops.

