



Item number 1999905376



Pants, ABENA Pants, Special M-L2, Premium

- ✓ Secure and comfortable fit with legs and extra elastics
- ✓ Super soft - for optimal comfort
- ✓ Support healthy skin balance
- ✓ Leakage barriers for extra protection
- ✓ Top-Dry Technology for dry surface to protect the skin
- ✓ 100% breathability for healthy skin
- ✓ Odour Control



Product Description

ABENA Pants Special is our Premium range of disposable pants that can be worn like normal underwear for moderate incontinence. The products have legs and unique elasticity in both the legs and crotch area providing the ultimate fit during physical activities.

The products features the Top-Dry Technology for optimum comfort and skin integrity. The modern pant is made from breathable material and are fitted with leakage barriers for extra protection. When combined, these functionalities keep the skin in balance and maintains a natural pH skin balance. The elastic threads offers a secure and comfortable fit and the odour control system minimises unpleasant smells.



Specifications

Base name	Pants
Brand	ABENA
Sub-brand	Pants
Size	Special M-L2
Color	White
Features	Light blue colourcode
Material	Breathable textile backsheet (BTBS), ECF fluff, elastane, nonwoven, PE, PP, resin, SAP
Length/depth	810 mm
Width	760 mm
Hip size	80-135 cm
Tolerance	+/- 10%
Absorption Rothwell	1700 ml
Quality level	Premium
Certifications	CE, MD, UKCA
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1, ISO 13485, ISO 14001, ISO 9001
Directives, regulations and acts	MDR (EU) 2017/745
Safety Instructions and Warnings	Keep out of reach of children.
Because of the products composition it has an expiry date from the production date on:	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Dispose with household waste.

Packaging Disposal Instructions

Can be recycled or incinerated.

Instructions for use/application

Disposable pant for moderate incontinence that can be worn as normal underwear. Put on the pants as you would underpants, without touching the inside surface of the pad. After use and when changing, tear open the sides. When both sides of the pants are torn open, roll it up and discard it. Product size recognition by colour codes, names and numbering system.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cII	120 pcs	78.6 cm	38.5 cm	37 cm	5703538458613
pck	18 pcs	26 cm	18 cm	35 cm	5703538458477



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



The product does not contain latex.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.



Odour control is a critical factor in ensuring discreteness for the user. ABENA's incontinence products are designed and engineered to ensure that minimal odour arises from the product.



Effective leakage barriers, both visible and invisible and latex-free elastics ensure a snug fit, which reduces the risk of leakages.



Unique Top-Dry Technology that provides a rapid absorption and ensures a dry surface even after multiple voids.



A combination of 100% breathability, Top-Dry and leakage protection provides a dry surface and healthy skin even after multiple voids maintaining the skin's natural pH balance.



Drop system indicating absorption. The more drops, the higher the absorption. Ranges from 1 to 9 drops.

