



Item number 1000021315

ABENA®

Shaped pad, ABENA San, 7/XXL, purple colourcode, Premium

- ✓ Longer pad to suit a range of body shapes
- ✓ High performance core with Top-Dry Technology
- ✓ Visible and invisible leakage barriers for extra protection
- ✓ Wetness indicator with graduation scale
- ✓ Support healthy skin balance



Product Description

ABENA San 7/XXL and 11/XXL are ABENA's range of anatomically shaped pads designed for bariatric use. The longer pad provides greater security to suit a range of body shapes. The high performance core with Top-Dry Technology ensures fast absorption and a dry surface even after multiple voids, maintaining skin integrity. The fully breathable pads are fitted with standing leakage barriers combined with built-in cross barriers at the rear. When combined, these functionalities keep the skin in balance and maintains a natural pH skin balance.

The odour control system minimises unpleasant smells. The wetness indicator with graduation scale guides the user and caretaker to the product with the right absorption level and product changes. The dermatologically tested products feature the Nordic Eco Label.

These bladder protection pads should always be used with special fixation underwear, like our Abri-Fix Pants Super 4-6XL products, to ensure a close fit and reduce the risk of leakage. Two-piece systems, a fixation pant and an incontinence pad, are a cost-effective solution, as they use less materials, have less impact on the environment and increase quality of life.



Specifications

Base name	Shaped pad
Brand	ABENA
Sub-brand	San
Size	7/XXL
Color	White
Features	Purple colourcode
Material	Breathable textile backsheet (BTBS), ECF fluff, elastane, nonwoven, PE, PP, resin, SAP
Length/depth	840 mm
Width	460 mm
Tolerance	+/- 10%
Absorption Rothwell	2000 ml
Gender	Unisex
Quality level	Premium
Certifications	CE, MD, Nordic Swan Ecolabel, Dermatologically tested, UKCA
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1, ISO 13485, ISO 14001, ISO 9001
Directives, regulations and acts	MDR (EU) 2017/745
Because of the products composition it has an expiry date from the production date on:	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.

Product Disposal Instructions	Dispose with household waste.
Packaging Disposal Instructions	Can be recycled or incinerated.

Instructions for use/application

Anatomically shaped pads for moderate to heavy incontinence specially designed for bariatric use due to the additional pad length. Easy product recognition by colour codes, names and numbering system.

These bariatric pads should always be used with special fixation underwear, like our Abri-Fix Pants Super 4-6XL. These help keep the pad in place, which ensures a close fit and reduces the risk of leakages. Pull fixation pants half-way up the legs. Fold the pad lengthways to create a bowl shape. Pass the pad between the legs from front to back. Position the upper rear end of the pad at the coccyx. Fan out the back of the pad horizontally to create a channel shape between the legs. Use the wetness indicator to check if the product is applied accurately in the center. Fan out the pad horizontally at the front. Pull fixation pants up into position. See our application guide brochure for detailed information.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cII	104 pcs	57.6 cm	39.6 cm	36.7 cm	5713571000731
pck	26 pcs	38 cm	18 cm	26 cm	5713571000724



Nordic Swan Ecolabel is an official Scandinavian environmental label. Products with the Nordic Swan Ecolabel meet stringent environmental requirements with regard to raw material, production, use, documentation, and safety.



Your guarantee that our products have been tested for ingredients which can potentially cause skin irritation. Tests conducted by an independent third-party laboratory.



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



The product does not contain latex.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.



The unique structure of our high performance core is designed for fast dispersion and encapsulation of urine, which results in a dry surface even after multiple voids.



Odour control is a critical factor in ensuring discreteness for the user. ABENA's incontinence products are designed and engineered to ensure that minimal odour arises from the product.



Effective leakage barriers, both visible and invisible and latex-free elastics ensure a snug fit, which reduces the risk of leakages.



Unique Top-Dry Technology that provides a rapid absorption and ensures a dry surface even after multiple voids.



Wetness indicator is a liquid-sensitive strip that changes colour or fades out indicating when it is time to change the product.



A combination of 100% breathability, Top-Dry and leakage protection provides a dry surface and healthy skin even after multiple voids maintaining the skin's natural pH balance.



Drop system indicating absorption. The more drops, the higher the absorption. Ranges from 1 to 9 drops.

