

Product Description

ABENA Pants Special is our Premium range of disposable pants that can be worn like normal underwear for moderate incontinence. The products have legs and unique elasticity in both the legs and crotch area providing the ultimate fit during physical activities.

The products features the Top-Dry Technology for optimum comfort and skin integrity. The modern pant is made from breathable material and are fitted with leakage barriers for extra protection. When combined, these functionalities keep the skin in balance and maintains a natural pH skin balance. The elastic threads offers a secure and comfortable fit and the odour control system minimises unpleasant smells.



Specifications

| _ | | | | |
|--|--|--|--|--|
| Base name | Pants | | | |
| Brand | ABENA | | | |
| Sub-brand | Pants | | | |
| Size | Special S-M2 | | | |
| Color | White | | | |
| Material | Breathable textile backsheet (BTBS), ECF fluff, elastane, nonwoven, PE, PP, resin, SAP | | | |
| Ingredients / Composition | Blue colourcode | | | |
| Length/depth | 790 mm | | | |
| Width | 760 mm | | | |
| Hip size | 60-110 cm | | | |
| Tolerance | +/- 10% | | | |
| Absorption Rothwell | 1700 ml | | | |
| Quality level | Premium | | | |
| Certifications | CE, MD, UKCA | | | |
| CE Class (Medical Devices) | Class I | | | |
| Product or test standards | ISO 11948-1, ISO 13485, ISO 14001, ISO 9001 | | | |
| Directives, regulations and acts | MDR (EU) 2017/745 | | | |
| Safety Instructions and Warnings | Keep out of reach of children. | | | |
| Because of the products composition it has an expiry date from the production date on: | 5 years | | | |
| Storage Instructions | Store dry, room temperature and no direct sunlight. | | | |
| Product Disposal Instructions | Dispose with household waste. | | | |
| | | | | |



Packaging Disposal Instructions

Can be recycled or incinerated.

Instructions for use/application

Disposable pant for moderate incontinence that can be worn as normal underwear. Put on the pants as you would underpants, without touching the inside surface of the pad. After use and when changing, tear open the sides. When both sides of the pants are torn open, roll it up and discard it. Product size recognition by colour codes, names and numbering system.

Packaging data

| Unit | Contains | Length | Width | Height | EAN |
|------|----------|---------|---------|--------|---------------|
| cll | 120 pcs | 78.6 cm | 38.5 cm | 37 cm | 5703538458460 |
| pck | 20 pcs | 26 cm | 18 cm | 35 cm | 5703538458422 |



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.











