



## TECHNICAL DATA SHEET

### **PRODUCT DESCRIPTION**

<b>PRODUCT CODE</b>	493CE03
<b>GENERAL DESCRIPTION</b>	MRSA KIT+1REGUL.+2UNIV.APPL.
<b>DETAILED DESCRIPTION</b>	10 Vipak each containing 50 film-film pouches each containing 2 flocked universal pink swabs and 1 flocked regular swab coated in a nutrient solution and 1 tube with screw wite cap and 1 ml of LQ Amies medium

### **PRODUCT COMPONENTS**

<b>MEDIUM</b>	<b>COMPOSITION</b>
LIQUID AMIES MODIFIED	Sodium Chloride Tween 80 Sodium Thioglycollate Sodio fosfato Calcium Chloride Beef extract Disodium Phosphate Distilled water

<b>SAMPLING TUBE</b>	<b>TUBE</b>	<b>DIMENSIONAL SPECIFICATIONS</b>
		TOTAL LENGTH: 81.30 mm ± 0.20 mm TUBE DIAMETER: 13.00 mm ± 0.20 mm INTERNAL DIAMETER: 11.50 mm ± 0.20 mm
<b>CAP</b>	<b>WHITE CAP</b>	<b>DIMENSIONAL SPECIFICATIONS</b>
		CAP LENGTH: 15.00 mm ± 0.20 mm CAP DIAMETER: 16.80 mm ± 0.20 mm

This bulletin (Rev.02-10/04/2018) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



## TECHNICAL DATA SHEET

<b>APPLICATOR</b>	UNIVERSAL PINK FLOCKED APPLICATOR	<b>DIMENSIONAL SPECIFICATIONS</b>
		FLOCKED TIP DIAMETER: 5.50 mm $\pm$ 0.50 mm TOTAL LENGTH: 151.0 mm $\pm$ 2.0 mm LENGTH OF THE PLASTIC SHAFT: 150.0 mm $\pm$ 1.0 mm PLASTIC TIP DIAMETER: 3.00 mm $\pm$ 0.10 mm STICK DIAMETER HANDLE PART: 2.50 mm $\pm$ 0.10 mm LENGTH OF THE FIBRE TIP: 16.0 mm $\pm$ 2.0 mm
<b>APPLICATOR</b>	FLOCKED REGULA SWAB	<b>DIMENSIONAL SPECIFICATIONS</b>
		STICK LENGTH: 150.0 mm $\pm$ 1.0 mm DIAMETER HANDLE PART: 2.50 mm $\pm$ 0.10 mm PLASTIC TIP DIAMETER: 3.00 mm $\pm$ 0.05 mm TOTAL LENGTH: 151.0 mm $\pm$ 2.0 mm BREAK POINT: 78.50 mm $\pm$ 1.00 mm FIBRE TIP LENGTH: 16.0 mm $\pm$ 2.0 mm FIBRE TIP DIAMETER: 5.5 mm $\pm$ 0.5 mm

### **PACKAGING SPECIFICATION**

<b>OUTER BOX QTY</b>	500
<b>PACKAGING</b>	FILM-FILM POUCH
<b>VI-PAK QUANTITY</b>	50

### **IMPORTANT INFO**

<b>STERILIZATION METHOD</b>	GAMMA IRRADIATION
<b>STORAGE TEMPERATURE</b>	5/25°C-DO NOT FREEZE
<b>SHELF LIFE</b>	15 MONTHS
<b>MEDIUM pH</b>	7.2 $\pm$ 0.3

This bulletin (Rev.02-10/04/2018) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.